

## SPECIFICATION SHEET

<b>Product Name:</b>	Mouse anti human IgG (H&L)McAb	<b>Catalog #:</b>	L1I00301
<b>Description:</b>	Mouse Antibody to human Immunoglobulin G (IgG), Heavy & Light chains		
<b>Specificity:</b>	Reacts with heavy chains on human IgG and with the light chains on all human immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal human Serum.		
<b>Host Animal:</b>	Mouse		
<b>Immunogen:</b>	Purified human IgG		
<b>Format:</b>	Affinity Purified, Liquid		
<b>Purification:</b>	>90% pure (SDS-PAGE). Immunoaffinity purified. Product is 0.2 micron filtered.		
<b>Concentration:</b>	5.3mg/ml	<b>Lot:</b>	L180601
<b>Buffer:</b>	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Specific methodologies have not been tested using this product.		
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

### FOR RESEARCH USE ONLY