

SPECIFICATION SHEET

Product Name:	Human AFP	Catalog #:	L2C002
Description:	Purified Human Alpha Fetoprotein (AFP) AFP is a glycoprotein containing about 4% carbohydrate. AFP is produced by the fetal liver and is not normally expressed in adults.		
Source:	Human fetal cord serum		
Format:	Purified, Liquid.		
Purification:	Column chromatography		
Concentration&Lot:	Activity: 1 mg/ml	Lot:	L160101
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4, 0.05%NaN3		
Preservative:	0.05%NaN3		
Applications:	Specific methodologies have not been tested using this product.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HbsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		