

SPECIFICATION SHEET

Description:	Recombinant Human GRP
Catalog #:	L2C03001
MW:	17kd
Source:	Recombinant, Pichia pastoris
Format:	Purified, Liquid
Purification:	Purified, > 90% (determined by SDS gelelectrophoresis)
Net Weight:	0.2mg
Lot:	L181101
Buffer:	Lyophilized from 10mM Tris-HCL, PH8.0 ± 0.2
Preservative:	None
Storage:	Store at -20 C. Avoid multiple freeze/thaw cycles.