

SPECIFICATION SHEET

Product Name:	Hepatitis B Surface Ag (HBs Ag)	Catalog #:	L1H00301
Description:	Hepatitis B Surface Antigen		
Source:	Human infected plasma		
Format:	Purified, Liquid		
Purification:	Ion exchange chromatography		
Concentration:	1.75 mg/ml	Lot. :	L181201
Buffer:	0.01M PBS, pH 7.4		
Preservative:	0.02% Sodium azide		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Inactivation:	Heated for 10 hours at 60°C before ultracentrifugation.		
Warnings:	All human source materials have tested negative for HIV1, HIV2 and HCV antibodies. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		
For Research Use Only.			