

SPECIFICATION SHEET

Product Name:	HBc Ag	Catalog #:	L1H00701
Description:	Recombinant Hepatitis B c Ag Hepatitis B c Antigen (Hepatitis B core Antigen, 144aa,14kD)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	>90% pure (SDS-PAGE) Ion exchange chromatography		
Concentration:	2.6 mg/ml	Lot.:	L160201
Buffer:	0.02M Sodium phosphate, pH 7.2		
Preservative:	none		
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

For Research Use Only