

SPECIFICATION SHEET

Product Name:	HIV-1 P24	Catalog #:	L2H002
Description:	HIV-1 P24 Recombinant antigen, Antigen Grade		
Source:	Recombinant ,E.coli		
Format:	Purified, Protein solution in Tris-HCl pH8.0,		
Purification:	Purity >90%		
Concentration&Lot:	2mg/ml	Lot:	L181201
Buffer:	20Mm Tris-HCl, pH 8.0 ± 0.02,		
Preservative:	None		
Application:	Immunoassays or Immunization.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		