

SPECIFICATION SHEET

Product Name:	McAb-HIV-1P24	Catalog #:	L2H00503
Description:	MAb to HIV-1 p24 Monoclonal Antibody to Human Immunodeficiency Virus Type 1 (HIV-1), p24 specific		
Specificity:	HIV-1 p24 (Recombinant antigen)		
Clone:	MLP3E12		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Recombinant Antigen		
Format:	Purified, Liquid		
Concentration:	1.8mg/ml (OD280nm)	Lot:	L190401
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	Suitable as a detection antibody in ELISA. Each laboratory should determine an optimum working titer for use in its particular application.		
Storage:	Store at -20°C. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.