

SPECIFICATION SHEET

Product Name:	Rabbit anti-cTnI PcAb	Catalog #:	L3C00403
Description:	Rabbit Anti Human cTnI Polyclonal Antibody,		
Specificity:	Recognizes an epitope on the cTnI Antigen.		
Clone:	None		
Host Animal:	Rabbit		
Isotype:	IgG		
Source:	serum		
Immunogen:	Human cTnI Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE).		
Concentration&Lot:	6.5 mg/ml	Lot :	L190501
Buffer:	10Mm PBS, pH 7.4 ± 0.2		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		