

SPECIFICATION SHEET

Product Name:	McAb to NT-Pro-BNP	Catalog #:	L3C00503				
Description:	Monoclonal Antibody to Human NT-pro-Brain Natriuretic Peptide (NT-proBNP)						
Specificity:	NT-proBNP						
Clone:	NB003						
Host Animal:	Mouse						
Isotype:	IgG						
Source:	Ascites						
Immunogen:	Recombinant NT-pro-Brain Natriuretic Peptide (NT-proBNP)						
Format:	Purified, Liquid						
Purification:	Protein G chromatography						
Concentration&Lot:	1mg/ml	Lot :	L160501				
Buffer:	10mM PBS, pH 7.4						
Preservative:	none						
Applications:	Can be used in an immunoassay for specific determination of human myoglobin. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L3C00503</td><td>L3C00504</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L3C00503	L3C00504
<u>Capture</u>	<u>Detection</u>						
L3C00503	L3C00504						
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.						