

SPECIFICATION SHEET

Product Name:	PcAb-sST2	Catalog #:	L3C01203
Description:	PcAb to sST2 antibody Rabbit Anti Soluble ST2 Polyclonal Antibody,		
Specificity:	Recognizes an epitope on the sST2 Antigen.		
Host Animal:	Rabbit		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Human sST2Antigen		
Format:	Affinity Purified, Liquid,Protein A chromatography		
Purification:	≥90% pure (SDS-PAGE).		
Concentration:	9mg/ml (OD280nm)	Lot :	L170901
Buffer:	0.01M PB, 150mm NaCl,pH7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	store at -20°C.		