

## SPECIFICATION SHEET

<b>Product Name:</b>	C-Peptide Antibody	<b>Catalog #:</b>	L3I00603
<b>Description:</b>	PcAb to C-Peptide antibody Rabbit Anti Human C-Peptide Antibody,		
<b>Specificity:</b>	Recognizes an epitope on the C-Peptide Antigen.		
<b>Clone:</b>	None		
<b>Host Animal:</b>	Rabbit		
<b>Isotype:</b>	IgG		
<b>Source:</b>	serum		
<b>Immunogen:</b>	Human C-Peptide Antigen		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	≥90% pure (SDS-PAGE). Protein A chromatography		
<b>Concentration &amp; Lot.:</b>	2.45mg/ml (OD280nm) (Lot. specific)		
<b>Buffer:</b>	PBS, pH 7.4 ± 0.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage &amp; Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		