

SPECIFICATION SHEET

Catalog #:	L1B00202	Lot #:	C140801				
Description:	MAb mouse anti-human S-100β Protein Monoclonal Antibody to S-100β Protein (S-100BB)						
Clone:	3A9						
Host Animal:	Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from Balb/c mice.						
Isotype:	IgG ₁						
Source:	Ascites						
Immunogen:	Purified human S-100 Protein						
Format:	Purified, Liquid						
Purification:	Protein A chromatography						
Concentration:	1.17 mg/ml						
Buffer:	PBS, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs in ELISA: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1B00201</td><td>L1B00202</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1B00201	L1B00202
<u>Capture</u>	<u>Detection</u>						
L1B00201	L1B00202						
Storage:	store at -20°C.						