

SPECIFICATION SHEET

| Catalog #: | L1B00203 | Lot #: | L180901 | | | | |
|-----------------------|--|---------------|---------|----------------|------------------|----------|----------|
| Description: | Rabbit anti-human S-100β Protein PcAb Rabbit anti human S-100β Protein (S-100BB) Polyclonal Antibody | | | | | | |
| Clone: | N/A | | | | | | |
| Host Animal: | Rabbit. | | | | | | |
| Source: | Serum | | | | | | |
| Immunogen: | Purified human S-100 Protein | | | | | | |
| Format: | Purified, Liquid | | | | | | |
| Purification: | Protein A chromatography | | | | | | |
| Concentration: | 4 mg/ml | | | | | | |
| Buffer: | PBS, pH 7.4 | | | | | | |
| Preservative: | None | | | | | | |
| Applications: | Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs in ELISA: <table><thead><tr><th><u>Capture</u></th><th><u>Detection</u></th></tr></thead><tbody><tr><td>L1B00201</td><td>L1B00203</td></tr></tbody></table> | | | <u>Capture</u> | <u>Detection</u> | L1B00201 | L1B00203 |
| <u>Capture</u> | <u>Detection</u> | | | | | | |
| L1B00201 | L1B00203 | | | | | | |
| Storage: | store at -20°C. | | | | | | |