

SPECIFICATION SHEET

<b>Product Name:</b>	Rabbit anti Recombinant AB1-42 Protein	<b>Catalog #:</b> L1B01003
<b>Description:</b>	PcAb to Recombinant AB1-42 Protein Rabbit anti-Amyloid beta peptide(AB1-42) Polyclonal antibody	
<b>Specificity:</b>	Reacts with Recombinant AB1-42 Protein.	
<b>Clone:</b>	N/A	
<b>Host Animal:</b>	Rabbit	
<b>Source:</b>	Serum	
<b>Immunogen:</b>	Recombinant AB1-42 Protein.	
<b>Format:</b>	Purified, Liquid	
<b>Purification:</b>	Protein A chromatography	
<b>Concentration&amp;Lot:</b>	3.2 mg/ml	<b>Lot:</b> L190301
<b>Buffer:</b>	10mM PBS, pH 7.4	
<b>Preservative:</b>	none	
<b>Applications:</b>	Suitable for use in immunoassay,	
<b>Storage&amp;Shelf life:</b>	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.