

SPECIFICATION SHEET

| | | | |
|-----------------------|---|---------------|---------|
| Catalog #: | L1V00103 | Lot #: | L180301 |
| Description: | Rabbit anti-25-OH-Vitamin D3-BSA (Calcifediol-BSA) Polyclonal Antibody | | |
| Specificity: | Recognizes an epitope on the 25-OH-Vitamin D3-BSA (Calcifediol-BSA) | | |
| Immunogen: | 25-OH-Vitamin D3-BSA (Calcifediol-BSA) | | |
| Host Animal: | Rabbit | | |
| Source: | Serum | | |
| Format: | Purified, Liquid | | |
| Concentration: | 3.3mg/ml | | |
| Buffer: | 0.01M PBS, pH 7.4 | | |
| Preservative: | None | | |
| Applications: | ELISA | | |
| Storage: | Store at -20 deg C for long term storage. Avoid repeated freeze-thaw cycles. | | |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.