

SPECIFICATION SHEET

Product Name: McAb to S100-A8 **Catalog #:** L6B00402

Description: Mouset anti Recombinant Human S100 Calcium Binding Protein A8

Specificity: Recognizes an epitope on Human S100 Calcium Binding Protein A8

Immunogen: Recombinant Human S100 Calcium Binding Protein A8

Host Animal: Mouse

Source: Serum

Format: Purified, Liquid

Concentration: 5.8mg/ml **LOT:** L180801

Buffer: 0.01M PBS, pH 7.4

Preservative: None

Storage&Shelf life: Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.