

SPECIFICATION SHEET

| | | | |
|--------------------------------|--|-------------------|---------|
| Product Name: | B19 VP1+2 | Catalog #: | L4H002 |
| Description: | Recombinant human B19 VP1+VP2 Antigen (B19 VP1+VP2) | | |
| Molecular Weight: | 33KD (calculated from sequence) 33KD (determined by SDS gelelectrophoresis) | | |
| Source: | E. coli | | |
| Format: | Purified, Liquid | | |
| Purification: | > 90% (determined by SDS page) | | |
| Concentration&Lot: | 2.5 mg/ml | Lot : | L170701 |
| Buffer: | 0.01M PBS pH 7.4 ± 0.2 0.5mM EDTA | | |
| Preservative: | None | | |
| Application: | Immunoassays as EIA, IFA or Immunization. | | |
| Storage&Shelf life: | Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles. | | |