

SPECIFICATION SHEET

Product Name:	Human anti HBcAg PcAb	Catalog #:	L1H00602
Description:	Human anti Hepatitis B c Antigen (HBcAg) Polyclonal Antibody		
Specificity:	HBcAg.		
Clone:	N/A		
Host Animal:	Human		
Isotype:	IgG		
Source:	Human Serum		
Immunogen:	HBcAg		
Format:	Purified, Liquid		
Purification:	>90% pure. Affinity purified.		
Concentration:	3.0mg/ml	Lot. :	L180301
Buffer:	10MM PBS, pH 7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.