

SPECIFICATION SHEET

Catalog #:	L3C00704	Lot #:	L190301
Description:	PcAb to D-dimer antibody Rabbit Anti Human D-dimer Antibody,		
Specificity:	Recognizes an epitope on the D-dimer Antigen.		
Clone:	None		
Host Animal:	Rabbit		
Isotype:	IgG		
Source:	serum		
Immunogen:	Human D-dimer Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE).		
Concentration:	6.5mg/ml (OD280nm)		
Buffer:	PBS, pH 7.4 ± 0.2		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	store at -20°C.		