

Description: L-T₄ L-Thyroxine Sodium Salt Pentahydrate

Catalog #: L2G01002

Molecular Formula: C₁₅H₁₀I₄NO₄Na·5H₂O

MW: 888.93

Purity: 80%

Lot: L160401

Forma: Powder

Storage: Store at 2 ~ 8°C , Against light or wet.

Dissolved: Dissolve one by one drops of NaOH 0.1mol/L, and make it dissolve completely. The solution is diluted to the appropriate concentration with 50% propylene glycol solution.(Note: this method is only applicable to the calibration of the product, for reference. Other uses of the dissolution method please refer to the relevant information.)

FOR RESEARCH USE ONLY.