

## SPECIFICATION SHEET

<b>Catalog #:</b>	LS1A01101	<b>Lot #:</b>	L160301				
<b>Description:</b>	McAb to APOE antibody Mouse anti Human APOE monoclonal Antibody,						
<b>Specificity:</b>	Recognizes Human APOE.						
<b>Host Animal:</b>	Mouse						
<b>Immunogen:</b>	Human APOE						
<b>Format:</b>	Liquid						
<b>Buffer:</b>	0.01M PBS, 0.15M NaCl, pH 7.4,						
<b>Concentration:</b>	5mg/ml						
<b>Applications:</b>	Suitable for use in IHC, IFA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><thead><tr><th><u>Capture</u></th><th><u>Detection</u></th></tr></thead><tbody><tr><td>LS1A01101</td><td>LS1A01102</td></tr></tbody></table>			<u>Capture</u>	<u>Detection</u>	LS1A01101	LS1A01102
<u>Capture</u>	<u>Detection</u>						
LS1A01101	LS1A01102						
<b>Storage:</b>	Store at -20°C. Avoid multiple freeze/thaw cycles.						

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.