

SPECIFICATION SHEET

Product Name:	Mouse anti human HAMA	Catalog #: L1I006
Description:	Mouse Antibody to human HAMA	
Specificity:	human HAMA	
Host Animal:	Mouse	
Immunogen:	Purified HAMA	
Format:	Liquid	
Purification:	>80% pure (SDS-PAGE).	
Concentration:	45mg/ml	Lot.: L190801
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2	
Preservative:	None	
Applications:	Specific methodologies have not been tested using this product.	
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.	

FOR RESEARCH USE ONLY